



# AJL OCULAR

## Perforated Conformer

<i>Brand name</i>	<b>AJL OCULAR</b>		
<i>References</i>	<b>CF-S-AJL</b>	<b>CF-M-AJL</b>	<b>CF-L-AJL</b>
<i>Size</i>	Small: 22 x 19	Medium: 24 x 21	Large: 25 x 22
<i>Description</i>	<p>The perforated conformer is a medical device used after surgical procedures such as enucleation or evisceration.</p> <p>The perforated conformer has a concave shape that fits the eye socket and has two holes to aid the drainage of secretions and the administration of drops and medication.</p>		
<i>Indications</i>	<p>The AJL perforated conformer is designed to maintain the conjunctival sac after enucleation or evisceration without hindering eyelid closure. The conformer maintains the separation between the tarsal and bulbar conjunctivae to prevent adhesion between the two tissues, while simultaneously preserving space and preventing cavity contraction.</p> <p>Its use is indicated after evisceration or enucleation surgery.</p>		
<i>Material</i>	It is made of a highly resistant and biocompatible thermoplastic translucent material of polyethylene and copolyester (PTEG).		
<i>Sterilization</i>	Ethylene oxide (ETO).		
<i>Supply</i>	Single-unit, sterile packaging.		
<i>Shelf life</i>	4.5 years.		

CE 0318

[ajlsa.com](http://ajlsa.com)

