



AJL OCULAR

Ocular Shield

<i>Brand name</i>	AJL OCULAR		
<i>References</i>	PROCUL-S-AJL	PROCUL-M-AJL	PROCUL-L-AJL
<i>Size</i>	Small: 22 x 19 x 1.7 mm	Medium: 24 x 21 x 1.7 mm	Large: 25 x 22 x 1.7 mm
<i>Description</i>	The AJL ocular shield is a medical device specially designed to protect the eyeball.		
<i>Indications</i>	It is designed to protect the globe during eyelid surgery.		
<i>Material</i>	High-strength PETG with a surface finish suitable for eye use.		
<i>Sterilization</i>	Ethylene oxide (ETO).		
<i>Supply</i>	Single-use, sterile packaging. Includes a latex-free suction cup for easy application and removal of the eye protector.		
<i>Shelf life</i>	4.5 years.		

CE 0318

ajlsa.com