



# AJL OCULAR

## Symblepharon Ring

<i>Brand name</i>	<b>AJL OCULAR</b>		
<i>References</i>	<b>SF-S-AJL</b>	<b>SF-M-AJL</b>	<b>SF-L-AJL</b>
<i>Size</i>	Small: 22 x 19	Medium: 24 x 21	Large: 25 x 22
<i>Description</i>	The AJL symblepharon ring is an oculoplastic device designed to prevent conjunctival fornix retraction. The central width of the hole prevents corneal contact while preserving the cornea.		
<i>Indications</i>	The symblepharon ring is designed to prevent retraction of the conjunctival fornix while preserving the cornea. It is placed under the patient's eyelid to prevent the development of a condition known as symblepharon (a partial or complete adhesion between the inner surface of the eyelid and the eyeball resulting from scarring). This condition can affect both eyelids, although it is more common in the lower eyelid.		
<i>Material</i>	It is made of a highly resistant, translucent thermoplastic material of polyethylene and copolyester (PTEG) and is biocompatible.		
<i>Sterilization</i>	Ethylene oxide (ETO).		
<i>Supply</i>	Single-unit, sterile packaging.		
<i>Shelf life</i>	4.5 years.		

CE 0318

[ajlsa.com](http://ajlsa.com)

